

## Laboratory Reference ranges

1. Open **Health One Configuration Tool**
2. **File**  
**Open**  
**Term**
3. Choose term (ie WBC) and highlight
4. Change **Entry Style** to Numeric
5. Change **Display Style** to Lab #1
6. Click **Apply**, and close window  
(Message: Do you want to save) **Yes**
7. **File**  
**Open**  
**Term**
8. Highlight your chosen term (ie. WBC) again
9. Click on **Default value** and enter **0**
10. A window will open. Enter the required reference range in the appropriate cell.
11. Click **OK**
12. Click **Apply** (The content has been changedí í í ..)
13. Click **Yes**
14. Re-enter Health one records and click **Tools** on the tool bar.
15. **Reload Dictionaries**
16. Click **OK**